



DESCRIÇÃO DOS CILINDROS AÇO



TABELA DE DIMENSÕES

Item	Modelos	Vol. Hidráulico	Cap. M3	Diâmetro	Altura	Peso
1	3.6Lts	3.6Litros	0.54	139.7	400mm	5
2	03Lts	03Litros	0.45	114.3	400mm	4.8
3	05Lts	05Litros	0.75	139.7	450mm	7.9
4	07Lts	07Litros	1	165.1	457mm	9.6
5	10Lts	10Litros	1.5	165.1	620mm	12.2
6	16Lts	16Litros	2.4	177.0	880mm	17.8
7	20Lts	20Litros	4	219.0	720mm	30.0
8	40Lts	40Litros	8	219.0	1320mm	48.0
9	50Lts	50Litros	10	235.0	1420mm	56.0

WWW.ROMED.COM.BR

ROMED EQUIPAMENTOS MÉDICOS
PABX: +55(11) 3718-1000
VENDAS@ROMED.COM.BR

